




## EPIC ClearView Product Requirements Trace Matrix - Device

**ORIGINATOR:**

Name	Job Title	Signature	Date
Scott Pletzer	Director of Quality & Program Management		03APR12

**APPROVALS:**

Name	Job Title	Signature	Date
Nancy Rizzo	Founder and CEO		4/3/12
Andrew Mason	Chief Technology Officer		4/3/12

**REVISION HISTORY**

Revision	Revision Date	Reason for Revision/Change Request	Revised By
000	03Apr12	Initial Release	S. Pletzer

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-1	Camera and lens assembly capable of capturing an energized image.	Lens assembly	Phillips Camera (SPC900NC/27) ClearView BOM, 06690 CCD Camera and Collar, 1282-1003 Westech Lens, 1282-1004	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Verification (PADT), #07612 Device Design Validation, ENG-017		
DR-2	Camera capable of image capture at a minimum of 15 pictures/10 sec exposure	Established by electrophotography techniques	Phillips Camera (SPC900NC/27) ClearView BOM, 06690	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Verification (PADT), #07612 Device Design Validation, ENG-017		
DR-3	Camera must be able to be controlled through the software for all available settings pertaining to image quality.	Established by electrophotography techniques	Phillips Camera (SPC900NC/27) ClearView BOM, 06690	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Verification (PADT), #07612 Device Design Validation, ENG-017 Camera Functions Validation Report, ENG-010 Application Functions Validation Report, ENG-009		
DR-4	A transparent electrode which can be configured for receiving a fingertip of a subject, such as on a top surface thereof. A transparent electrode able to withstand high voltages in excess of 25k volts at high frequencies. Able to generate a localized electromagnetic field around the finger.	Established by electrophotography techniques	Electrode Assembly, 06617 PCBA, Control Board, 06596	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Verification (PADT), #07612 Device Design Validation, ENG-017 Capture Functions Validation Report, ENG-011		

## EPIC ClearView Product Requirements Trace Matrix - Device

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-5	IEC 60601-1 compliant	Regulatory requirements	ClearView BOM, 06690	IEC Certification  Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016	IEC 60601-1 Report, ENG-002		
DR-6	IEC 60601-1-2 compliant				Device Design Verification (PADT), #07612		
DR-7	FCC, Part 15 compliant				Device Design Validation, ENG-017 Receiving Inspection Record, QAF-012		
DR-8	Device is designed to ensure proper grounding of electrical current to prevent patient and operator exposure to electrical currents.	Safety features identified			EPIC ClearView IEC 60601-1-2:2007 Final Test Report, ENG-002 EPIC ClearView FCC Rule Part(s) 15, Subpart B Final Test Report, ENG-003 EPIC ClearView FCC Rule Part(s) 15.20, 15.209 Final Test Report, ENG-004		
DR-9	Per IEC 60601-1, device materials, manufacture and use do not present a flammability risk.	Flammability risks		Device Design Verification and Validation			
DR-10	Device operation is repeatable and reliable against IFU stated conditions.	User preference	ClearView BOM, 06690	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Verification (PADT), #07612 Device Design Validation, ENG-017 DV-002 Clinical Study Report		
DR-11	Device has an EPIC branded logo on the top surface.	EPIC marketing requirement	Housing Cover, Finished, 06646	Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016	Receiving Inspection Records, QAF-012		



## EPIC ClearView Product Requirements Trace Matrix - Device

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-12	Device is usable by patients regardless of variability in size and position (e.g., ambidexterity, prone, etc.) and moveable to accommodate patient access without special modifications.	Patient anthropometric requirements and patient comfort	ClearView BOM, 06690	Documentation System Policy, QA-005 Device Design Verification and Validation Clinical Trial Performance Monitoring	Device Design Validation, ENG-017 EPIC ClearView IFU, CS-100 EPIC ClearView System Performance Reports: • GBMC Clinical Trial, ENG-018 • DV-001 Clinical Trial, ENG-019 • DV-002 Clinical Trial, ENG-020 • EPIC-003 Clinical Trial, ENG-021		
DR-13	Device can be shown to capture images that accurately represents the subject's galvanic skin response	Reliability requirements	ClearView BOM, 06690	Documentation System Policy, QA-005 Device Design Verification and Validation	Device Design Validation, ENG-017 EPIC ClearView IFU, CS-100 EPIC ClearView Clinical Trial Reports: • GBMC Clinical Trial Report • EPIC ClearView Algorithm Development Report • EPIC-003 Clinical Trial Report		
DR-14	PER IPX0 rating, device must resist the ingress of liquids, and be operated under standard room conditions, similar to other solid state electronic devices.	Device use, care, and storage requirements	ClearView BOM, 06690	IEC Certification Clinical Trial Performance Monitoring	IEC 60601-1 Report, ENG-002 EPIC ClearView IFU, CS-100 EPIC ClearView System Performance Reports: • GBMC Clinical Trial, ENG-018 • DV-001 Clinical Trial, ENG-019 • DV-002 Clinical Trial, ENG-020 • EPIC-003 Clinical Trial, ENG-021		

## EPIC ClearView Product Requirements Trace Matrix - Device

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-15	Device/system compatible with IEC 60601-1, IEC 60601-2, and FCC Part 15 requirements	User expectation	ClearView BOM, 06690	IEC Certification Clinical Trial Performance Monitoring	Final Test Reports: <ul style="list-style-type: none"><li>• EPIC ClearView IEC 60601-1-2:2007 Final Test Report, ENG-002</li><li>• EPIC ClearView FCC Rule Part(s) 15, Subpart B Final Test Report, ENG-003</li><li>• EPIC ClearView FCC Rule Part(s) 15.20, 15.209 Final Test Report, ENG-004</li></ul> EPIC ClearView System Performance Reports: <ul style="list-style-type: none"><li>• GBMC Clinical Trial, ENG-018</li><li>• DV-001 Clinical Trial, ENG-019</li><li>• DV-002 Clinical Trial, ENG-020</li><li>• EPIC-003 Clinical Trial, ENG-021</li></ul>		
DR-16	Maintain certifications and regulatory requirements: 21 CFR 801.415 - Maximum Acceptable Level of Ozone. Device cannot be released for use without these in place.	Applicable safety standards	ClearView BOM, 06690	Ozone emissions testing Risk Analysis Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016	EPIC ClearView Ozone Emissions Verification, ENG-005 EPIC Scanner Risk Analysis Documentation Form, DHD-101 EPIC ClearView Failure Modes and Effects Analysis, DHD-102 Receiving Inspection Record, QAF-012		
DR-17	Device materials must be non-toxic wherever skin contacts the device.	Toxicity risks	ClearView BOM, 06690	Inspection and Release of Incoming,	Electrode Assembly, 06617		

## EPIC ClearView

### Product Requirements Trace Matrix - Device

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-18	Device materials are biocompatible with user and/or patient.	Biocompatibility		In-process and Finished Devices and Other Materials, QA-016  Regulatory Review	Housing Base, Painted, 06647 Housing Cover, Finished, 06646 FDA 510k Submission		
DR-19	Device glass electrode surface needs to be cleanable with 70% IPA and not damaged between patient uses.	Product life and maintenance requirements	Cleaning Cloth, 1282-1500 70% Isopropyl Alcohol Swabs, 1282-1504	Device Design Verification and Validation Documentation System Policy, QA-005 Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016	Device Design Verification (PADT, #07612) Device Design Validation, ENG-017 EPIC ClearView IFU, CS-100 Receiving Inspection Record, QAF-012 Calibration Testing Reports, ENG-006 & ENG-007 EPIC ClearView System Performance Reports: <ul style="list-style-type: none"><li>GBMC Clinical Trial, ENG-018</li><li>DV-001 Clinical Trial, ENG-019</li><li>DV-002 Clinical Trial, ENG-020</li><li>EPIC-003 Clinical Trial, ENG-021</li></ul>		
DR-20	Device needs to be calibratable per instructions for proper performance.			Calibration Probe, 06680 Calibration Probe Cable, 06683 EPIC ClearView IFU, CS-100	Clinical Trial Performance Monitoring		
DR-21	Device has a power on switch accessible by the user which clearly indicates the device is powered on.	Access control	ClearView Device, Packaged BOM, EPC-100 Power Entry Module	Documentation System Policy, QA-005 Device Design	Device Design Verification, XXX-XXX Device Design Validation, XXX-XXX		



## EPIC ClearView

### Product Requirements Trace Matrix - Device

Device (DR)							
Ref. #	Requirement	Source of Requirement	Design Inputs (Specification)	Verification Methodology	Design Outputs (Verification Results)	Validation Results	Status
DR-22	Device housing must be secured in such a way that tampering will not occur without tools.		Assembly, 06686  Rear Panel LED Cable, 06670  EPIC ClearView System Description, SD-003  EPIC ClearView Software Description, SD-001  USB Hardware Key, 1282-1515  EPIC ClearView IFU, CS-100	Verification and Validation  Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016	EPIC ClearView IFU, CS-100  Receiving Inspection Record, QAF-012		
			Documentation System Policy, QA-005  Inspection and Release of Incoming, In-process and Finished Devices and Other Materials, QA-016  EPIC ClearView Software Validation	Validation Reports: <ul style="list-style-type: none"><li>• Camera Functions Validation Report, ENG-010</li><li>• Capture Functions Validation Report, ENG-011</li><li>• Patient Demographics Validation Report, ENG-012</li><li>• Users Functions Validation Report, ENG013</li><li>• Report Functions Validation Report, ENG-014</li></ul>			
DR-23	Device capable of communicating with a computer.						
DR-24	Device complies with applicable international and national standards and requirements.	<ul style="list-style-type: none"><li>• 21 CFR Part 820</li><li>• 21 CFR Part 11</li><li>• ISO 13485</li><li>• MDD</li></ul>	EPIC Quality System	Documentation System Policy, QA-005  Internal and External Audits	Internal Audit Reports, QAF-008  FDA 510k Submission  ISO 13485/MDD External Audit		